



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Dental X-Ray Unit

**Product Code :**  
SAMCOMEDICAL1400001A006



**Description :**

Dental X-Ray Unit ½ The dental X-ray package is a true insurance of quality with all the benefits of the Advanced Technocracy Inc ½ An extremely high dielectric coefficient is obtained through the use of specially treated oil ½ The high voltage transformer comes with the best insulators ½ A special expansion chamber compensates for the dilation of the oil caused by heat from the tube ½ The protection against leakage radiation is stricter than the strictest international standards. The hooded tube with grid, compensates for voltage fluctuation ½ Input voltage: 230V AC single phase Timer: Electric with digital display of mAs Output: 8 mA - 70 KVP Stand: Floor model/Wall model-ray tube: 1.5mm focal Spot Cone: Parallel Dental X-Ray Unit, Dental Instrument exporters, Dental Instrument manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),