



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Dental Compressor Oil Less

**Product Code :**  
SAMCOMEDICAL1400001A007



**Description :**

Dental Compressor Oil Less 1/2 HP Dental compressor seal unit fitted with safety valve for air release excess sorted, auto cut off system pressure gauge and drain valve air rotor air compressor mounted on high pressure tank entirely powder coated paint Power: 3/4HP Capacity: 80L/min Speed: 1400rpm Noise: 58dB(A)/40" Tank: 38L Pressure: 116psi Weight: 28kgs Dimension: 42x42x62cm Packing size: 47x47x70cm Guarantee: 24 month Dental Compressor Oil Less, Dental Instrument exporters, Dental Instrument manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),