



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ultrasonic Scaler

Product Code :
SAMCOMEDICAL1400001A010



Description :

Ultrasonic Scaler Scaling Function: Ultrasonic Scaler is a powerful piezoelectric scaler with stable ultrasonic oscillation of 32KHz with very wide energy scope Titanium Hand piece is autoclavable Great patient comfort with high frequency and low amplitude Handpiece: never heat up after long time use Efficient performance to remove intractable calculus Fine tuning control on power supply and water spray Endodontic Function: Turn the power below half of TOP SELECTOR and replace tip with endo kit (Endoholder = U-file) to operate unit for root canal preparation. Smoother canal walls grant better obturation Rapid and superior operating for root canal treatment Root canal enlargement and swirling irrigation for debridement Minimal dentin removal keeps structural integrity of the tooth Technical Specifications: Power Supply 230V~ or 115V~ (on request) 50/60 Hz Frequency 32 KHz 1/2 Output Power to Tip Operator adjustable from 20 to 40 watts Weight 2.2Kg Dimension W150 X D130 X H110mm Water Pressure 1~4Kg / cm2 Fuse F0.5A (for 230V unit) / FIA (for 115V unit) Standard Accessories Optional Accessories 1 Handpiece & Cable Code Item 3 Scaling Tips COM315 Ti-Handpiece 1 Torque Wrench COM316-A1 Ti No.1 (Normal curve) 1 Pedal COM316-A2 Ti No.2 (Normal curve) 1 Power Cord COM316-A3 Ti No.3 (Flat) Tips stand SE201 U-File#25 Aluminum box COM203-2 Torque Wrench Ultrasonic Scaler, Dental Instrument exporters, Dental Instrument manufacturers



Samco Medical

Website: www.samcomedical.com,