



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Curing Light

Product Code :
SAMCOMEDICAL1400001A011



Description :

Curing Light ½ Power source: 230V~115V or 115V~(on request), 50/60Hz ½ Weight: 2.4 Kgs. ½ Dimension: W200 x D140 x H80(mm) Curing Light Function: ½ Complete curing restoration with a 75 watt, 12V, halogen lamp ½ Light guide with 360 degree rotation and can be autoclavable ½ 3 seconds interval to turn on or off the lamp controlled by computer to protect and lengthen use life of lamp ½ Individual curing time selection with beep tone signal ½ Automatical turning off fan after cool down period ½ Provides fuse and thermal cut-off switch to prevent overheating Curing Light, Dental Instrument exporters, Dental Instrument manufacturers



Samco Medical

Website: www.samcomedical.com,