



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dental Double Walled Acrylizer

Product Code :
SAMCOMEDICAL1400001A014



Description :

Dental Double Walled Acrylizer is a double walled, inside made of stainless steel and out side mild steel sheet painted with attractive stove enamel. Temperature is controlled by hydraulic thermostat from 5°C above ambient to 100°C with a sensitivity of + 5°C. Acrylizer consists of Mechanical Timer can be set 1/2 Hour to 2 Hours, brass drain plug, control knob, On/Off switch. Acrylizer Inner Size : 300 X 250 X 350 mm. Dental Double Walled Acrylizer, Dental Instrument exporters, Dental Instrument manufacturers.



Samco Medical

Website: www.samcomedical.com,