



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Labour Table Foot End Folding

Product Code :
SAMCOMEDICAL9002A001

Description :

Labour Table Foot End Folding We are counted among the reputed organizations, engaged in offering our clients with Labour Table Foot End Folding. All our products are available in different sizes, shapes and designs at market leading prices to fulfill the demands and needs of customers. In addition to this, we stringently check all the products on various parameters to ensure flawlessness. Other Information: $\frac{1}{2}$ Overall Size: L 1830 mm x W 870 mm x Ht. 800 mm $\frac{1}{2}$ Frame: 50 mm x 25 mm $\frac{1}{2}$ Two section Stainless steel top with head low/up position with screw mechanism $\frac{1}{2}$ Detachable and folding foot end $\frac{1}{2}$ Height Adjustment and head low/up position by screw and lever mechanism $\frac{1}{2}$ $\frac{1}{2}$ U' notch with drainage tray $\frac{1}{2}$ Available in full stainless steel top /stainless steel top & epoxy Powder coated base Labour Table Foot End Folding, Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,