



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Labour Table Telescopic

Product Code :
SAMCOMEDICAL9002A002



Description :

Labour Table Telescopic ½ Overall Size L 1830 mm x W 870 mm x H 800 mm Frame 50 mm x 25 mm
Rectangular Stainless steel tube ½ Three section Stainless steel top ½ Back rest on ratchet; head low/up by
screw mechanism ½ Middle Section with 'U' notch with bucket ½ Sliding foot end ½ Railings on three sides
½ Available in full stainless steel top /stainless steel top & epoxy Powder coated base Labour Table
Telescopic, Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,