



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Mayo Instrument Tables

Product Code :
SAMCOMEDICAL9002A006



Description :

Mayo Instrument Tables The Mayo Instrument Table comes with a stainless steel tray and a tubular MS framework, mounted on four swivel castors. Fabricated from high grade raw materials, the Mayo Instrument Table can be availed at competitive market prices. Mayo Instrument Tables, Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,