



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hydraulic Operation Table (Electric)

Product Code :
SAMCOMEDICAL9002A007



Description :

Hydraulic Operation Table (Electric) *Electric operation table with x-y top longitudinal displacement movement finger touch remote control with radiolucent top * Electric Hydraulic Operation Table Specially designed for use with C-arm Image Intensifier. The entire uppertable top can be removed to entire side i.e to headend or to foot end side so that the body may be screened when the table base is in the standard position. eccentrically positioned table column and interchangeable head and leg section. Base and column cover made of stainless steel . complete with set of accessories including mattress. provide smooth, easy and accurate positioning controlled by noiseless actuators Remote control system with cable connection offers Electric Hydraulic Operation Table X-Y Top Longitudinal displacement movement * Electric Hydraulic Operation Table Up and down Position * Electric Hydraulic Operation Table Tendlenberg & Reverse Tendlenberg POSITION Electric Hydraulic Operation Table Lateral position * Electric Hydraulic Operation Table headside raising & lowering position Hydraulic Operation Table (Electric), Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,