



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hydraulic Operation Table (Electric)

Product Code :
SAMCOMEDICAL9002A007



Description :

Hydraulic Operation Table (Electric) *Electric operation table with x-y top longitudinal displacement movement finger touch remote control with radiolucent top * Electric Hydraulic Operation Table Specially designed for use with C-arm Image Intensifier. The entire upper table top can be removed to either side i.e. to headend or to foot end side so that the body may be screened when the table base is in the standard position. eccentrically positioned table column and interchangeable head and leg section. Base and column cover made of stainless steel . complete with set of accessories including mattress. provide smooth, easy and accurate positioning controlled by noiseless actuators Remote control system with cable connection offers Electric Hydraulic Operation Table X-Y Top Longitudinal displacement movement * Electric Hydraulic Operation Table Up and down Position * Electric Hydraulic Operation Table Trendelenburg & Reverse Trendelenburg POSITION Electric Hydraulic Operation Table Lateral position * Electric Hydraulic Operation Table headside raising & lowering position Hydraulic Operation Table (Electric), Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,