



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Hydraulic Operation Table (Mechanical)

**Product Code :**  
SAMCOMEDICAL9002A008



**Description :**

Hydraulic Operation Table (Mechanical) \* : Hydraulic Operation Table are Tubular framework of MS mounted on PVC stumps.  $\frac{1}{2}$ : Hydraulic Operation Table Stainless steel top in three sections  $\frac{1}{2}$  : Hydraulic Operation Table Trendelenburg/ Rev. Trendelenburg position adjustable by pneumatic gas spring system.  $\frac{1}{2}$ : Hydraulic Operation Table are Backrest section adjustable on ratchet.  $\frac{1}{2}$ : Hydraulic Operation Table are Leg end section can slide under the main section.  $\frac{1}{2}$  : Hydraulic Operation Table Height adjustable with a pair of knee crutches  $\frac{1}{2}$ : Hydraulic Operation Table Provided with side railings, hand grip & lithotomy rod. : Hydraulic Operation Table Size: 72"(L) x 27"(W) x 30"(H) (Approx) Finish: Pre treated to prevent rust & Epoxy powder coated Hydraulic Operation Table (Mechanical), Hospital Table exporters, Hospital Table manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),