



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hydraulic Operation Table (Mechanical)

Product Code :
SAMCOMEDICAL9002A008



Description :

Hydraulic Operation Table (Mechanical) * : Hydraulic Operation Table are Tubular framework of MS mounted on PVC stumps. $\frac{1}{2}$: Hydraulic Operation Table Stainless steel top in three sections $\frac{1}{2}$: Hydraulic Operation Table Trendelenburg/ Rev. Trendelenburg position adjustable by pneumatic gas spring system. $\frac{1}{2}$: Hydraulic Operation Table are Backrest section adjustable on ratchet. $\frac{1}{2}$: Hydraulic Operation Table are Leg end section can slide under the main section. $\frac{1}{2}$: Hydraulic Operation Table Height adjustable with a pair of knee crutches $\frac{1}{2}$: Hydraulic Operation Table Provided with side railings, hand grip & lithotomy rod. : Hydraulic Operation Table Size: 72"(L) x 27"(W) x 30"(H) (Approx) Finish: Pre treated to prevent rust & Epoxy powder coated Hydraulic Operation Table (Mechanical), Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,