



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hydraulic Obstetric Labour Table

Product Code :
SAMCOMEDICAL9002A010



Description :

Hydraulic Obstetric Labour Table * Hydraulic Obstetric Labour Table Framework made MS tubes. 1/2 Hydraulic Obstetric Labour Table ARE Legs fitted with PVC stumps. 1/2 Hydraulic Obstetric Labour Table MS sheet top in two sections. 1/2 Hydraulic Obstetric Labour Table Backrest section adjustable on ratchet. Hydraulic Obstetric Labour Table Size: 183(L) x 51(W) x 81(H) cms. (Approx) Finish: Pre-treated to prevent rust and Epoxy Powder Coated. Hydraulic Obstetric Labour Table, Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,