



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Examination Couch with Cabinet

Product Code :
SAMCOMEDICAL9002A011



Description :

Examination Couch with Cabinet *: Examination Couch with Cabinet Made from MS sheet, machine pressed.
1/2: Examination Couch with Cabinet Easy and instantaneous tilt of backrest section with the help of Gas spring system. 1/2: Examination Couch with Cabinet Foam padded, upholstered top for patient comfort. 1/2: Examination Couch with Cabine Cabinet in lower section for storage of medicines, bandages etc., using sliding doors for maximum utility within limited space. 1/2 : Examination Couch with Cabinet Three drawers in upper section and swinging tray for B.P. apparatus. : Examination Couch with Cabinet Size : 187(L) x 51(W) x 81(H) cms. (Approx) Backrest : 45L x 51W cms Main Frame : 138L x 51W cms Examination Couch with Cabinet, Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,