



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Gynaceo/Delivery Table

Product Code :
SAMCOMEDICAL9002A012



Description :

Gynaceo/Delivery Table Three sections with independent control and Trelenburg position. Padded and upholstered with washable and flame-retardant plastic material. Complete with stainless steel bowl, 16 1/2" 32 cm., and two adjustable anatomic knee crutches. Gynaceo/Delivery Table, Hospital Table exporters, Hospital Table manufacturers



Samco Medical

Website: www.samcomedical.com,