



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Gynae Examination Couch

Product Code :
SAMCOMEDICAL9003A001

Description :

Gynae Examination Couch Backed by our efficient workforce, we are able to provide our customers with Gynae Examination Couch. All our products are widely acclaimed by the customers for their corrosion & abrasion resistance, sturdy construction, low maintenance and longer service life. Further, we provide tailor-made solution for all the products as per the specifications & demands of customers. Other Information: Overall Size: L 1830 mm x W 575 mm x H 850 mm Frame: CRCA sheets and rectangular tubes 3 Cabinets; 1 with lock, drawer with lock Coir Mattress with 'U' notch Examination Trolley with three trays Pre-treated and epoxy powder coated Gynae Examination Couch, Hospital Couch exporters, Hospital Couch manufacturers



Samco Medical

Website: www.samcomedical.com,