



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Examination Couch

Product Code :
SAMCOMEDICAL9003A002



Description :

Examination Couch We bring forth for our customers a wide range of Examination Couch that are available in different sizes and models at market leading prices. All the products are designed in accordance with international quality standards, using premium raw material. Furthermore, we thoroughly check all the products on various parameters to ensure flawlessness. Other Information: Overall Size: L 1830 mm x W 575 mm x H 850 mm Frame: CRCA Sheets and Rectangular tubes 100 mm thick form mattress 4 Cabinets; 1 with lock Drawer with lock Sliding stepping tool & writing pad pre-treated and epoxy powder coated Examination Couch, Hospital Couch exporters, Hospital Couch manufacturers



Samco Medical

Website: www.samcomedical.com,