



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Deluxe Bedside Locker

Product Code :
SAMCOMEDICAL9005A001

Description :

Deluxe Bedside Locker Overall size: 40L x 40W x 82H cms Complete locker box made of 22G CRC machine pressed sheet closed from three sides and a suitable drawer for extra storage Stainless steel top of 22G in matt finish Complete locker mounted on two PVC Stumps in front and two 5 cms. dia casters at the back for easy moving Pre-treated and epoxy powder coated finish Deluxe Bedside Locker, Hospital Lockers exporters, Hospital Lockers manufacturers



Samco Medical

Website: www.samcomedical.com,