



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Deluxe Bed Side Locker

**Product Code :**  
SAMCOMEDICAL9005A006



**Description :**

Deluxe Bed Side Locker We are engaged in offering our clients with Deluxe Bed Side Locker that are available in different sizes, shapes and designs at market leading prices. All the products are widely acclaimed by the customers for their corrosion & abrasion resistance, dimensional accuracy and robust construction. Moreover, we have a team of expert quality controllers, which checks all the products on certain parameters to ensure flawlessness. Other information: Overall size: W 370 mm x D 305 mm x H 800 mm Two stainless steel tops with 3 sides raised and front turned down Railings on 3 sides above main cabinet Handle and magnetic latch on door CRCA steel box mounted on 20 mm x 20 mm CRCA sq. tube frames Pre-treated and epoxy powder coated Deluxe Bed Side Locker, Hospital Lockers exporters, Hospital Lockers manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),