



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Plastic Sprayed Bed Side Lockers

Product Code :
SAMCOMEDICAL9005A009



Description :

Plastic Sprayed Bed Side Lockers The Hospital Plastic Sprayed Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Plastic Sprayed Bedside Locker has a stainless steel top and mounted on two swivel castors. Hospital ABS Bedside Locker. The Hospital ABS Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Bedside Locker has a stainless steel top and mounted on two swivel castors. Plastic Sprayed Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers



Samco Medical

Website: www.samcomedical.com,