



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Economy Model Folding Wheel Chair

**Product Code :**  
SAMCOMEDICAL9006A002



**Description :**

Economy Model Folding Wheel Chair 4" Caster Powder Coated Seat and back easily removable and replaceable by any household person Fine and durable upholstery for seat and back Plastic Hand grips Padded arm rest Most Comfortable Heavy duty construction 3/4" Tube Mild Steel Powder Coated Solid tyre wheels Rear wheel locks Metal Side Panel Aluminium leg rest Economy Model Folding Wheel Chair, Folding Wheel Chair exporters, Folding Wheel Chair manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),