



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electronic Wheelchair

Product Code :
SAMCOMEDICAL9006A011



Description :

Electronic Wheelchair & Wheelchair Electronic With Safety Lamp & Adjustable Armrest Pad, Detachable Sideboard & Adjustable Angle Backrest & Maximum Distance covered on Single Charge: 20 km & Speed Range: 1-6 km/hour & Climbing Power: 12 & Battery: 28AH, 12V x 2, Motor: 168W x 2 Electronic Wheelchair, Folding Wheel Chair exporters, Folding Wheel Chair manufacturers



Samco Medical

Website: www.samcomedical.com,