



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hydraulic Operation Table

Product Code :
SAMCOMEDICAL9007A001

Description :

Hydraulic Operation Table We are instrumental in offering our clients with Hydraulic Operation Table that are mainly used for all major and minor operating procedures e.g. general surgery, gynaecology - obstetrics, urology, cardiothoracic and ophthalmic. These products are manufactured using high-grade raw material, which is sourced from the reliable vendors of the market. Conform to international quality norms & standards, our products are known for their corrosion & abrasion resistance, sturdy construction, low maintenance and zero defect. It is supplied with the following standard accessories: 1/2 S. S. Arm Rest 1 No 1/2 Anaesthetic Screen 1 No. 1/2 Lithotomy Leg Holders with stirrups 1 Set 1/2 Leather Wristlets 1 Set 1/2 Padded Leg Rest (Gutter Type) can be supplied at an extra cost The table will attain the following positions: 1/2 Trendelenburg 1/2 Reverse Trendelenburg 1/2 Lateral on either sides 1/2 Flex 1/2 Reflex 1/2 Chair It will have the following salient features: 1/2 Detachable head rest is easily adjustable to any desired position, above or below table top 1/2 Detachable leg section facilitates ortho-attachment to be assembled at later stage. It is also most convenient to gynaecologist & urologist 1/2 Table top can be rotated 360 1/2 through base 1/2 Kidney bridge is provided to attain 13 cms height 1/2 All corners of base and table top are rounded to avoid any injury Hydraulic Operation Table, Operation Theater Tables exporters, Operation Theater Tables manufacturers



Samco Medical

Website: www.samcomedical.com,