



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Operating Table

**Product Code :**  
SAMCOMEDICAL9007A002



**Description :**

Operating Table i.e. Versatile Operation Table for general surgical procedures i.e. Accurate positions can be achieved with easy mechanical operation i.e. Five Section Radiolucent Top for X-ray Imaging i.e. Heavy Duty Double Cylinder Hydraulic Pump i.e. Base and Column covered with Stainless Steel fittings for Hygiene and Durability i.e. Movable on High Quality Casters i.e. Equipped with Base Locking Mechanism for firm rigid locking of the table i.e. Interchangeable Leg and Head Section Standard Accessories: i.e. 01 Shoulder Supports (1 pair) i.e. 02 Lateral Supports ( 1 pair) i.e. 03 Knee Crutches/Leg Holder Straps ( 1 pair) i.e. 05 Arm Board ( With radiolucent top) ( 1 pair) i.e. 07 Anaesthetic Screen ( 1 piece) Technical Specification: i.e. Length: 1860 mm i.e. Width: 510 mm i.e. Maximum Height: 980 mm i.e. Minimum Height 760 mm i.e. Trendelenberg: Upto 30 i.e. Reverse Trendelenberg: Upto 30 i.e. Lateral Tilt: 20 i.e. Both side i.e. Head Rest: 20 i.e. Up 90 i.e. Down Operating Table, Operation Theater Tables exporters, Operation Theater Tables manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),