



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Orthopedic Extension Device For Operating Table

Product Code :
SAMCOMEDICAL9007A003



Description :

Orthopedic Extension Device For Operating Table This orthopedic Extension Device can be used with any Operating Table as it independently stands on the floor
• Specially Designed for Intra Operative X-Ray Photos and Image Intensifiers use
• Removable Pelvis Rest and X-Ray Translucent Pelvis Stem
• Stirring Instruments allow Flexion/Raising of Legs and Stirring Ball Joint Locked from Foot End
• Two Traction Devices On sliding Spot with Ball Joints Allowing Flexions and Rotations of Foot with Traction Shoe
• Padded orthopedic knee Crutches provides support to the legs of the patient
• Requires Small Storage space as it is easily collapsible
Orthopedic Extension Device For Operating Table, Operation Theater Tables exporters, Operation Theater Tables manufacturers



Samco Medical

Website: www.samcomedical.com,