



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
4 Reflector Light with Stand

Product Code :
SAMCOMEDICAL9008A001

Description :

4 Reflector Light with Stand ÷ Dome Head :515mm Dia ÷ Halogen Bulbs :4nos. 12v x 50 watt ÷ Light intensity at 1 mt. :70,000 Lux ÷ Intensity Control :4 Stage ÷ Height Adjustment :600 mm approx ÷ Action Radius :1250mm ÷ Possible Movements :Radial, Angular & Axial ÷ Color Temperature :4200 K ÷ 300K ÷ Temp.rise in field :3÷6 c from Amb.Temp ÷ Low voltage unit with CVT :1 No. (C.V.T optional) ÷ Input Supply Voltage :220 volt A.C ÷ 10% 50Hz ÷ 2% 4 Reflector Light with Stand, Operation Theater Lights exporters, Operation Theater Lights manufacturers



Samco Medical

Website: www.samcomedical.com,