



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Operation White Light 7 Reflector

Product Code :
SAMCOMEDICAL9008A004



Description :

Operation White Light 7 Reflector We have gained expertise in providing our customers with Operation White Light 7 Reflector. These products are widely acclaimed by the customers for their corrosion & abrasion resistance, sturdy construction, low maintenance and excellent performance. In addition to this, we pack all the products using premium packaging material to ensure complete safety during transportation. Dimensions: Dome Head:715mm Dia. Halogen Bulbs:7nos. 12v x 50 watt Light intensity at 1 mt.:1,20,000 Lux Intensity Control:4 Stage Height Adjustment:500mm Action Radius:1850mm Possible Movements:Radial Angular & Axial Color Temperature:4200 K \pm 1/2 300K Temp rise in field:3 \pm 1/2 - 6 \pm 1/2 c from Amb. Temp Low voltage unit with CVT:1 No. (C.V.T optional) Input Supply Voltage:220voltA.C \pm 1/2 10% 50Hz \pm 1/2 2% Operation White Light 7 Reflector, Operation Theater Lights exporters, Operation Theater Lights manufacturers



Samco Medical

Website: www.samcomedical.com,