



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Operation Lights

Product Code :
SAMCOMEDICAL9008A005



Description :

Operation Lights A) Dome Head :715mm Dia. $\frac{1}{2}$ Halogen Bulbs: 7nos. 12v x 50 watt $\frac{1}{2}$ Light intensity at 1 mt.: 1,20,000 $\frac{1}{2}$ Lux Intensity Control: 4 Stage $\frac{1}{2}$ Height Adjustment: 500mm Action $\frac{1}{2}$ Radius: 1850mm $\frac{1}{2}$ Possible Movements: Radial, Angular & Axial Color Temperature :4200 K $\frac{1}{2}$ 300K $\frac{1}{2}$ Temp.rise in field: 4 $\frac{1}{2}$ - 8 $\frac{1}{2}$ c from Amb.Temp $\frac{1}{2}$ Low voltage unit with CVT: 1 No. (C.V.T optional) $\frac{1}{2}$ Input Supply Voltage: 220volt A.C $\frac{1}{2}$ 10% 50Hz $\frac{1}{2}$ 2% B) Dome Head :515mm Dia. $\frac{1}{2}$ Halogen Bulbs: 4nos. 12v x 50 watt $\frac{1}{2}$ Light intensity at 1 mt.: 70,000 $\frac{1}{2}$ Lux Intensity Control: 4 Stage $\frac{1}{2}$ Height Adjustment: 550mm Action $\frac{1}{2}$ Radius: 1550mm $\frac{1}{2}$ Possible Movements: Radial, Angular & Axial Color Temperature :4200 K $\frac{1}{2}$ 300K Temp.rise in field :3 $\frac{1}{2}$ -6 $\frac{1}{2}$ c from Amb.Temp $\frac{1}{2}$ Low voltage unit with CVT:1 No. (C.V.T optional) $\frac{1}{2}$ Input Supply Voltage: 220volt A.C $\frac{1}{2}$ 10% Operation Lights, Operation Theater Lights exporters, Operation Theater Lights manufacturers



Samco Medical

Website: www.samcomedical.com,