



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
S. S. Stool

Product Code :
SAMCOMEDICAL9016A002



Description :

S. S. Stool We are counted among the reputed organizations, engaged in offering our clients with S. S. Stool. All our products are available in different sizes, shapes and designs at market leading prices to fulfill the demands and needs of customers. In addition to this, we stringently check all the products on various parameters to ensure flawlessness. Highlights: $\frac{1}{2}$ Framework made of Stainless Steel $\frac{1}{2}$ Top surface is made of Stainless Steel $\frac{1}{2}$ Size: 320 mm Dia Top x 370 mm Dia Base x 540mm height Type: $\frac{1}{2}$ With rexine top $\frac{1}{2}$ Without rexine top S. S. Stool, Hospital Bedside Revolving Stools exporters, Hospital Bedside Revolving Stools manufacturers



Samco Medical

Website: www.samcomedical.com,