



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Operating Lights

Product Code :
SAMCOMEDICAL4008A003



Description :

Ceiling Operating Lights Ceiling Operating Lights are 19½ dia metallic dome fitted with halogen bulb (1x 24V/150W) & heat absorbing colour correcting filters Ceiling Operating Lights are Provided with plastic diffusers for shadow reduction. Spring balanced Ceiling Operating Lights are Provided with sterilizable handle and low voltage supply unit Ceiling Operating Lights, LED OT Light exporters, LED OT Light manufacturers



Samco Medical

Website: www.samcomedical.com,